Certificate of Analysis

Freeze-Dried Reference Serum – Target titer ranges*

LOT

RS10

14 May 2021

May 2026

AE ELISA	2500 – 4800
ART ELISA	600 – 2400
IBD ELISA	4000 – 10000
IBV ELISA	1900 – 5100
NDV ELISA	3800 – 9200
REO ELISA	2500 – 6800

Note:

**Titer range data has been formulated from multiple runs with different kit lots.
Please ensure this batch is included in your version of the BioChek software.*

Freeze-Dried Reference Serum

For use as a control serum in
BioChek Poultry Antibody ELISAs

Instructions for use

REF CR300

LOT RS10

BioChek (UK) Ltd
Unit 5 Kings Ride Park
Kings Ride, Ascot, Berkshire SL5 8BP
United Kingdom
E-mail: info@biochek.com
Website: www.biochek.com

Product description and intended use

The reference serum can be used for monitoring day-to-day variation in poultry antibody ELISAs and is produced by BioChek (UK) Ltd. The reference serum contains significant and known antibody levels against several major avian pathogens. A Certificate of Analysis is provided, giving either target titers or titer ranges for the different assays used.

Contents

- Each vial contains 0.20 ml (200 µl) serum in a freeze-dried presentation. Store the vials at 2-8 °C / 36-46 °F until use.
- Instructions for use (this leaflet), including Certificate of Analysis with target titers.

Precautions

- The control serum has been heat treated for 30 minutes at 56 °C / 133 °F to reduce the risk of possible (heat labile) infectious components. For extra safety handle all biological materials as potentially infectious.
- Treat all biological materials as potentially biohazardous. Wear gloves, safety glasses and a laboratory coat. Decontaminate used plates and waste including washings with a strong oxidising agent before disposal. Refer to product Safety Data Sheet for additional information.
- Use calibrated micropipettes and clean or sterile pipette tips for making all serum dilutions.
- Do not eat, drink or smoke in the laboratory. Do not pipette by mouth.
- Do not use control serum past its expiration date and do not mix different lot numbers.

Reconstitution of Freeze-Dried Reference Serum

- Use only one vial of Freeze-Dried Reference Serum at a time.
- Reconstitute the reference serum with 0.20 ml (200 µl) of BioChek green diluent, using a calibrated micropipette and a clean or sterile pipette tip.
- Allow to stand for approximately 15-30 minutes and fully dissolve the contents by gentle swirling of the vial. After complete dissolution of the contents mix by vortexing the vial to produce a homogenous serum.
- Reconstituted reference serum (when kept at 2-8 °C / 36-46 °F and stored in an upright position) can be used for up to 8 weeks for the testing of dilution techniques. After 8 weeks remains of undiluted reference serum should be discarded.

Using the Freeze-Dried Reference Serum

- Thoroughly mix the contents of the reference serum by pipetting up and down with a clean pipette tip or by vortexing, but try to avoid foaming as it will make accurate pipetting difficult.
- Determine the number of plates to be tested and set up the appropriate number of dilution plates. If samples are to be tested in duplicate, also use duplicates for the reference serum.
- Dilute the reference serum in exactly the same way and at the same time as all the other samples you are testing.
- Dispense the required amount of diluted reference serum in selected wells. At least 2 wells per plate should be used.


Evaluating and monitoring results

- The control serum titer value obtained can be compared with the control serum titer value stated on the Certificate of Analysis for each BioChek poultry ELISA kit.
- When the obtained titer value is within the expected ranges, one can be confident about quantitative comparability of results with other assays run in the same laboratory and other laboratories.
- Record test results of the control serum over time as titer groups. For more detailed information, endpoint titers may be used.
- Plot test results graphically to quickly monitor results over time, e.g. using BioChek software, or spreadsheet software or graphical paper.
- This Certificate of Analysis for Freeze-Dried Reference Serum only provides target titer ranges.

PLEASE NOTE: Target titers listed are only to be used as an indication of expected titer levels for each assay. Each laboratory has to establish its own target titer levels and monitor those over time.

Please ensure this batch is included in your version of the BioChek software.

For more information, please contact BioChek at:

BioChek B.V.
Fokkerstraat 14
2811 ER Reeuwijk
The Netherlands

Tel: +31 182 582 592
Fax: +31 182 599 360
e-mail: info@biochek.com
website: www.biochek.com

Symbol Glossary

Product code

Batch number


Manufacturing date

Manufacturer

Expiry date

Authorised representative in the European Community

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