

Certificate of Analysis

Ornithobacterium Rhinotracheale Antibody Test Kit

Product Code: CK108

Lot No: FS9465

Expiry Date: 25th October 2024

| | Kit Component | Lot No | | Kit Component | Lot No |
|---|-------------------|--------|---|------------------------|--------|
| 1 | Coated Plates | CH7212 | 6 | Sample Diluent Reagent | CH7899 |
| 2 | Conjugate Reagent | CH8158 | 7 | Wash Buffer Powder | CH6318 |
| 3 | Substrate Tablets | CH7952 | 8 | Negative Control | CH6861 |
| 4 | Substrate Reagent | CH8061 | 9 | Positive Control | CH8159 |
| 5 | Stop Solution | CH8053 | | | |

Quality Control Results

Mean Negative Control

Optical Density = 0.140

Mean Positive Control

Optical Density = 0.828

Pre-diluted Reference Control
(Product Code CR100, Lot no. RF18)

Antibody titre = 3718
(range: 2000 – 6500)

Date of test: 25th October 2023

The above result was generated by performing the test according to the package insert.

The above information is correct and forms part of the BioChek Quality Control for ORT ELISA kit product code CK108.



Released by QC – 27th October 2023



Andrew Russell
Quality Control Supervisor